


# Transanal Ileus Decompression Tube Set




Colonoscopic Retrograde Transanal Large Intestine  
Decompression & Drainage Tube Sets

A long nasotracheal tube is said to give best choice for conservative treatment of ileus. However, an emergency operation is often required for the left side of large intestine with cancerous ileus where decompression is not successful.

Transanally inserted, CLINY TRANSANAL ILEUS TUBE SET can effectively reduce pressure, which gives promising various advantages. This allows detailed preoperative examinations together with the conservative treatment suitable for disease to which emergency operation is not applicable.



Decompression tube



Guide wire



Dilator designed for opened forceps (DOF)



Dilator for dilating stenosis (DDS)

## Applications

Decompression and drainage of enteral content on the onset of cancerous ileus (left side of large intestine).

## Features and Advantages

### Facilitating transanal approach

Dedicated (specific) dilator for dilating stenosis (DDS) extends the stenosis into which decompression tube can be easily inserted.

### Producing higher decompression effect

Opened tip (end opening) and 4 side openings facilitate efficient suction of content. During suction, the air vent prevents excessive depressing pressure from acting on the intestine.

### Two-channel connector for easier irrigation and drainage

The two-channel connector allows easier and more hygienic irrigation and drainage procedure after a decompression tube is placed.

### Contrast rings for locating the balloon

The contrast rings at the both ends of the balloon are useful for locating the balloon.


### One-way valve for preventing the backflow of intestinal contents

The one-way valve prevents the backflow of intestinal contents through the air vent lumen.




Two-channel connector

## Basic Procedures




Insert the colonoscope immediately before reaching the stenosis of the large intestine.

Note: Be sure to follow the procedure under open endoscopy as well as X-ray fluoroscopy, determining the course of the product tip with care.




Insert the guide wire into the opening of the colonoscope for forceps to be passed through the stenosis.

Note: Never force the guide wire to be inserted because it can make perforation in the intestine.




Insert the dilator designed for the opened forceps (DOF) along the guide wire via the opening of the colonoscope for forceps to be passed through the stenosis.

Note: Be careful to operate the dilator avoiding risk of perforating the intestine with its tip.




Remove the colonoscope while the guide wire and DOF are left. Infuse olive oil into the inner surface of the DDS which is subsequently inserted along the DOF.

Note: Be careful to operate the dilator avoiding risk of perforating the intestine with its tip.



Remove the DDS while the guide wire and DOF are left. Fill the decompression tube with olive oil through its suction port.

Continue to insert the decompression tube along the DOF until the balloon part exceeds the stenosis.



Remove the guide wire and DOF and then infuse 30ml of sterilized water into the balloon.



## Specifications

### Transanal Ileus Decompression Tube Set

Product No.	Contents of the Set	
002 380 0220	①Transanal decompression tube	1pc 22Fr (7.3mm) L-1,200mm, Open tip, 4 side holes, Balloon capacity: 30mL
	②Guide wire	1pc 0.052" (1.32mm) L-3,000mm
	③Dilator designed for opened forceps (DOF)	1pc 8Fr (2.7mm) L-2,000mm
	④Dilator for dilating stenosis (DDS)	1pc 26Fr (8.6mm) L-1,000mm, Radiopaque tip type

With two-channel connector

Packaging: 1set per box, Sterilized with EtO Gas.

### Transanal Decompression Tube

002 385 0220	22Fr (7.3mm) L-1,200mm, Open tip, 4 side holes, Balloon capacity: 30mL
--------------	--

With two-channel connector

Packaging: 1pc per box, Sterilized with EtO Gas.

**Read the instructions for use thoroughly before use.**

- Please note that a partial amendment in specifications is made occasionally, as studies and developments have been performed continuously.



**CREATE MEDIC CO.,LTD.**

011

#### Head Office

8F, Shin-Yokohama Center Building, 2-5-15 Shin Yokohama, Kohoku-ku, Yokohama, Kanagawa, 222-0033 Japan

#### R & D Center - Export Division

3-25-6 Tonomachi, Kawasaki-ku, Kawasaki, Kanagawa, 210-0821 Japan

Tel: +81 44-577-0191 Fax: +81 44-577-0195

CE 0123