



MADE OF SILICONE

With excellent biocompatibility, it minimizes irritation to the urethral mucosa and eliminates concerns about latex allergies.

LARGE INNER LUMEN

The large inner lumen ensures sufficient flow rate.



EXCELLENT INSERTABILITY AND FLEXIBILITY

The balloon is highly flexible, making it gentle on the bladder, while the tube has an appropriate level of firmness, ensuring excellent insertability.

ALL SILICONE ANTIBACTERIAL FOLEY CATHETER

FEATURE

This is an antibacterial Foley Catheter with an antibacterial agent (silver) incorporated into the catheter tube. The results of tests for antibacterial activity and efficacy shows an antibacterial activity value of 2.0 or higher against Escherichia coli, Pseudomonas aeruginosa, MRSA, and Staphylococcus aureus, and its antibacterial effectiveness against these bacteria is confirmed.

Product Code	Size	Balloon Capacity	Valve Color
800 001 7145	12Fr	10mL	White
800 001 7146	14Fr		Green
800 001 7147	16Fr		Orange
800 001 7148	18Fr		Red
800 001 7149	20Fr		Yellow
800 001 7150	22Fr		Violet
800 001 7151	24Fr		Blue

Catheter size: Total length 430 mm Package: 10 pcs/box Sterilized