





SELF-CATHETERIZATION SET



Catheter
for
Urethral Self-Catheterization

CREATE MEDIC CO., LTD.



1. Foldable catheter case

Since the central part of the catheter case is foldable, it can be folded into a small size.

2. Hooks applicable to multiple purposes

A hook provided at the case upper end is designed to act as a fixing device when holding the case and as a hook to place the case on a door knob or similar devices while the catheter is in use.

3. Special cap keeps cleanliness

A cap with special structure does not only seal the case interior from exterior, but also allows sterilizer to freely enter the catheter interior, thereby always keeping the catheter clean.

4. Provided with an external container

An external container, provided as a standard accessory, is convenient for transporting the catheter case, and protects patient's privacy since the container conceals what is contained.

Self-Catheterization Set

Product No.	Catheter type	Size & Length	Side holes	Content
04-4110	Female	10Fr, L:165mm	1 hole	Catheter×1 Catheter case with hook ×1 Cap×1
04-4112		12Fr, L:165mm		
04-4114		14Fr, L:165mm		
04-4210	Male	10Fr, L:300mm		
04-4212		12Fr, L:300mm		
04-4214		14Fr, L:300mm		
04-4312	Male-L	12Fr, L:350mm		
800 000 5376		14Fr, L:385mm		
04-4409	Pediatric	9Fr, L:265mm	1 hole	
04-4512	Tieman	12Fr, L:300mm	1 hole,	
800 000 5375		14Fr, L:300mm	w/ Tieman Tip	
04-4614	C.U.R	14Fr, L:395mm	5 holes, with numeral markings	
04-4616		16Fr, L:395mm		
04-4618		18Fr, L:395mm		
04-4620		20Fr, L:395mm		

All Silicone Self Catheter (for exchange)

Product No.	Type	Size
04-5110	Female	10Fr
04-5112		12Fr
04-5114		14Fr
04-5210	Male	10Fr
04-5212		12Fr
04-5214		14Fr
04-5312	Male-L	12Fr
800 000 5478		14Fr
04-5409	Pediatric	9Fr
04-5512	Tieman	12Fr
800 000 5479		14Fr
04-5614	C.U.R	14Fr
04-5616		16Fr
04-5618		18Fr
04-5620		20Fr

Package: NOT sterilized, 5 sets in a box C. U. R= Continent Urinary Reservoir



External Case

Type	Size
Female Type	60×117×15
Male Type	60×190×15
Pediatric Type	
Tieman Type	
Male-L Type	60×232×15
C.U.R Type	



CREATE MEDIC CO.,LTD.

004

Head Office

8F, Shin-Yokohama Center Building, 2-5-15 Shin Yokohama, Kohoku-ku, Yokohama, Kanagawa, 222-0033 Japan

R & D Center - Export Division

3-25-6 Tonomachi, Kawasaki-ku, Kawasaki, Kanagawa, 210-0821 Japan

Tel: +81 44-577-0191 Fax: +81 44-577-0195

