


# Tierrett Catheter for cats (with stabilizer)

#### Urinary catheter for cats

(Registration No. of Medical Equipment for Animal Use: 28 Douyaku1210)

High reputation catheter for cats newly added a stabilizer. It is suitable for flushing, urethral catheterization/placement.



"Tierrett"

"Tierrett", derived from the German word

"Tierrettung" (aid for an animal), is a brand of CREATE MEDIC for animal

medical products that were developed with the aim of helping sick animals.







# Tierrett Features of catheter for cats (with stabilizer)



### Soft and flat stabilizer

The soft stabilizer, the shaft side is flat and gently fits to the cat's body



Soft and flat stabilizer

eature 2

## Two side holes closed tip

In addition to the conventional open tip type for flushing, a new lineup of closed tip type



Left:Open tip

Right: Closed tip

Feature 3

Choose 3 length types depending on cat's individual differences

In addition to the conventional 120 mm, the lineup of 80 mm and 100 mm



12typs of lineup (Refer to size composition)

### Product Specifications

Size composition

Size Compos Size designation	Tip shape	Effective length	Product cord	JAN cord	Note
3Fr with stabilizer	Open tip type	80mm	800-001-4619	4542187202908	Package : 5-pack box/ sterilized
		100mm	800-001-4620	4542187202915	
		120mm	800-001-4621	4542187202922	
Tube O. D.: 1.10mm	Closed tip	80mm	800-001-4622	4542187202939	2 side holes, with stylet Package : 5-pack box/ sterilized
		100mm	800-001-4623	4542187202946	
		120mm	800-001-4624	4542187202953	
4Fr with stabilizer	Open tip type	80mm	800-001-4625	4542187202960	Package : 5-pack box/ sterilized
		100mm	800-001-4626	4542187202977	
		120mm	800-001-4627	4542187202984	
Tube O. D.: 1.33mm	Closed tip type	80mm	800-001-4628	4542187202991	2 side holes, with stylet Package : 5-pack box/ sterilized
		100mm	800-001-4629	4542187203004	
		120mm	800-001-4630	4542187203011	



**CREATE MEDIC CO., LTD.** 

https://www.createmedic.co.jp/english/

[Contact information]

Animal Medical Business Department E-mail: oem-infoch@createmedic.co.jp

[Head Office] 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan



Control No: 28-1210-20170912