

Tierrett

Extension Tube

Tierrett 1 Port Extension Tube

Registration No. of Medical Equipment for Animal Use: 1 Douyaku 970






Tierrett Y type 2 Port Extension Tube

Registration No. of Medical Equipment for Animal Use: 30 Douyaku 3058







Product features

The Y type 2 port has a long tube and can be operated or mixed injection at a position away from the infusion site.

It has the same functions and physical properties as an extension tube for human, its safety is guaranteed.

Achieves a low price with animal-only products.

You can choose from 3 types.

Lineup

Product name	Tierrett 1 Port Extension Tube	Tierrett Y type 2 Port Extension Tube	
Port part image			
Type	T type connector	Straight connector	L type connector
Connector image			
Tube effective length	90mm	300mm	300mm
Product cord	800-001-6471	800-001-5989	800-001-5990
JAN cord	4542187221367	4542187216523	4542187216530

50-pack box/EOG sterilized

☐ Please note that the product specifications are subject to change without notice to reflect the results of our consistent R&D efforts.
☐ Read the package insert before using the product.



CREATE MEDIC CO., LTD.

https://www.createmedic.co.jp/english/

(Contact information)

Animal Medical Business Department E-mail: oem-infoch@createmedic.co.jp

(Head office) 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan



Control No: 30-3058-20200115