

Tierrett Guidewire

Registration No. of Medical Equipment for Animal Use: 28 Douyaku 2112



Offered special design for animals

"hydrophilic guidewire" at cheaper prices

Developed "for urinatoral stent placement"



OFeatures

- Two types are available choosing according to a case
- · Comes with hydrophilic coating
- Achieves both insertability and removability
 Comes with hydrophilic Ni-Ti cable core

"Tierrettung" (aid for an animal), is a brand of CREATE MEDIC for animal medical products that were developed with the aim of helping sick animals.

Tierrett

Guidewire



F e a

0

S

I

Two types are available: 0.016 inch (supported 2.1Fr stent), 0.018 inch (supported 2.5Fr stent)

I

High insertability is ensured by a moderately elastic hyper-elastic Ni-Ti (nickel titanium alloy) cable core and hydrophilic coating.

Ш

Adopted angle type to reduce tissue damage.

Product Specifications

Size designation	Outer diameter	Nose shape	Effective length	Product cord	JAN cord	Supported stent
0.016"	0.41mm	Angle	600mm	800-001-4718	4542187203899	Tierrett Ureteral Stent 2.1Fr
0.018"	0.46mm	Angle	600mm	800-001-4495	4542187201666	Tierrett Ureteral Stent 2.5Fr

1-pack box/EOG sterilized

 \square Please note that the product specifications are subject to change without notice to reflect the results of our consistent R&D efforts.

□Read the package insert before using the product.

Combination produc

Tierrett Ureteral Stent (Registration No. of Medical Equipment for Animal Use: 28 Douyaku 2112)

Size designation	Effective length								
2.1F r	80mm	100mm		120mm	140mm		160mm		
2.5F r 120mm		П	140mm			160mm			

1-pack box/EOG sterilized



CREATE MEDIC CO., LTD.

https://www.createmedic.co.jp/english/

[Contact information]

Animal Medical Business Department E-mail: oem-infoch@createmedic.co.jp

[Head Office] 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan

