Tierrett tracheal Tube for Large Animals

Registration No. of Medical Equipment for Animal Use: 2 Douyaku 2685



Desinged for foals and calves. Can be operated tracheal intubation with sufficient tube length. Also fit to use for large dogs.

Lineup

Size	Total length	I. D.	0. D.	Product cord	JAN cord
ID7.0	L-600mm	7.0mm	10.0mm(30Fr)	800-001-7053	4542187227192
ID9.0		9.0mm	12.0mm(36Fr)	800-001-7054	4542187227208
ID11.0		11.0mm	15.0mm(45Fr)	800-001-7055	4542187227215
ID13.0	L-650mm	13.0mm	18.0mm(54Fr)	800-001-7056	4542187227222
ID15.0	L-750mm	15.0mm	20.0mm(60Fr)	800-001-7057	4542187227239

1-pack box/Unsterilized

🎇 Please note that the product specifications are subject to change without notice reflecting the results of our consistent R&D efforts Read the package insert before using the product.

CREATE MEDIC CO., LTD.

A flexible body tube that made of silicone



CREATE MEDIC CO., LTD.

https://www.createmedic.co.jp/english/

(Contact information) Animal Medical Business Department E-mail : oem-infoch@createmedic.co.jp

(Head office) 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan



Control No: 2-2685-20210517