
## Designed exclusively for animals in pursuit of easy to use



## Tierrett Foley Catheter

Registration No. of Medical Equipment for Animal Use: 26 Douyaku 3980



## Recommend to large male dogs


For dogs whose urine leaks at 8Frl


## Lineup

mint.	A STATE OF THE STA	AND CHEST AND COURT	MARKET BARRETTE THE STATE OF	AD DIVERSION OF SHAPE OF		AND AN ARCHEOLOGICAL PROPERTY OF THE PROPERTY
	Size	O.D.	Balloon capacity	Effective length	Product cord	JAN cord
New	6Fr	2.0mm	1mL	650mm	800-001-3467	4542187191356
				850mm	800-001-3469	4542187191370
	8Fr	2.7mm	3mL	650mm	800-001-3470	4542187191387
				850mm	800-001-3472	4542187191400
	10Fr	3.3mm	3mL	650mm	800-001-6679	4542187223446
				850mm	800-001-6680	4542187223453

5-pack box/EOG sterilized

Please note that the product specifications are subject to change without notice to reflect the results of our consistent R&D efforts. Read the package insert before using the product.







Moderately elastic hydrophilic stylet for improve insertability



Easy-to-handle lightweight /compact joint



Reduced weight to 50% of conventional products and 55% of length





Multi-cap that can be used in 3 ways Catheter tip





Stopcock



CREATE MEDIC CO., LTD. https://www.createmedic.co.jp/english/

(Contact information)

Animal Medical Business Department E-mail: oem-infoch@createmedic.co.jp

[Head office] 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan



Control No: 30-1826-20201014