

Tierrett PEG Kit

Registration No. of Medical Equipment for Animal Use: 29 Douyaku 316








Features

O Kit Components

- 1)PEG bumper catheter
- 2Loop wire
- 3 Exclusive fixing plate
- **4**Connection 2 port adapter
- **5Clamp**
- @Drape
- 7Introducer cannula

with inner needle(14 G)

- **Syringe 10 mL**
- 9Gauze
- (1)Small incision scalpel
- 11) Disinfection sponge



OBumper type

The risk of accident 'self' removal is reduced.

OTransparent tube

You can easily check the feeding residue of lumen.

OExclusive fixing plate Since catheter is fixed with the exclusive fixing plate,

it is difficult to fall off and not in the way.



Product Specification

Size designation	Tube O.D.	Tube I.D.	Effective length	Bumper diameter	Product cord	JAN cord
16F r	φ5.3mm	φ3.0mm	700mm	21.5mm	800-001-5555	4542187212150
20F r	φ6.7mm	φ4.2mm	700mm	25mm	800-001-5028	4542187206807

1-Kit/box • EOG sterilized

- \square Please note that the product specifications are subject to change without notice to reflect the results of our consistent R&D efforts.
- □Read the package insert before using the product.



CREATE MEDIC CO., LTD.

https://www.createmedic.co.jp/english/

[Contact information]

Animal Medical Business Department E-mail: oem-infoch@createmedic.co.jp [Head Office] 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan



Control No:29-316-20181015