


Tierrett Ureteral Stent

Ureteral stent tube for cats

Registration No. of Medical Equipment for Animal Use: 28 Douyaku 2112

"Ureteral stent tube" for cats

Made in Japan "full lineup"



"Tierrettung" (aid for an animal), is a brand of CREATE MEDIC for animal medical products that were developed with the aim of helping sick animals.





Tierrett Ureteral Stent

O The tube is made of pliable
and elastic polyurethane

• High insertability
by delicate tip works

O High insertability is ensured by hydrophilic coating stent and dilator tube

Feature 2.

Feature 3.

O Eight types are available choosing according to a case from

both 2.1 Fr, 2.5 Fr tube size and length

Product Specifications

Size designation	O. D. I. D.	Accessories	Effective length	Product cord	JAN cord	Recommended guidewire (optional)
			80mm	800-001-4487	4542187201581	
0.15	O. D.: 0.71 I. D.: 0.45 Unit: mm	Dilater tube O. D.: 0.73mm	100mm	800-001-4488	4542187201598	Tierrett Guidewire O.O16"
2. 1Fr Green			120mm	800-001-4489	4542187201604	
Green			140mm	800-001-4490	4542187201611	
			160mm	800-001-4491	4542187201628	
0.55	O. D.: 0.83	D"	120mm	800-001-4492	4542187201635	Tierrett Guidewie
2. 5Fr Blue	I. D.: 0.53	Dilater tube O. D.: 0.86mm	140mm	800-001-4493	4542187201642	0.018"
Dide	Unit:mm	O. D 0.86mm	160mm	800-001-4494	4542187201659	0.018

1-pack box/EOG sterilized

□Please use together with our recommended guidewire "Tierrett guidewire" (optional).

□Please note that the product specifications are subject to change without notice to reflect the results of our consistent R&D efforts.

□Read the package insert before using the product.



https://www.createmedic.co.jp/english/

[Contact information]

Animal Medical Business Department E-mail: oem-infoch@createmedic.co.jp

[Head Office] 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan

