

Tierrett

Foley Catheter

Registration No. of Medical Equipment for Animal Use: 28 Douyaku 1077

Pursued “special design for animals” and
“a veterinarian usability”



○Features

- Comes with multi-cap that can be used in 3 ways
- Made of human-body-friendly silicon
- Lightweight compact joint

“Tierrett”, derived from the German word “Tierrettung” (aid for an animal), is a brand of CREATE MEDIC for animal medical products that were developed with the aim of helping sick animals.



Tierrett Foley Catheter

Features

- Useful **multi-cap** has **both syringe connection and plugging functions.**
 - It can be used in 3 ways as a catheter tip, lure tip, and plug.



- Easy-to-handle **lightweight/compact joint**



- Reduced weight to 50% of conventional company products and 55% of length

- Achieves insertability
Comes with **moderately elastic stylet**



- Made of human-body-friendly **silicon**
 - Hypoallergenic to urethral mucosa, reducing the adherence of urine constituent

Product Specification

Size designation	Outer diameter	Balloon capacity	Effective length	Product cord	JAN cord
6 F r	2.0mm	1mL	300mm	800-001-4611	4542187202823
8 F r	2.7mm	3mL	300mm	800-001-4612	4542187202830

5-pack box/EOG sterilized

- Please note that the product specifications are subject to change without notice to reflect the results
- Read the package insert before using the product.



CREATE MEDIC CO., LTD.

<https://www.createmedic.co.jp/english/>

【Contact information】

Animal Medical Business Department E-mail: oem-infoch@createmedic.co.jp

【Head Office】 2-5-25 Chigasaki Minami, Tsuzuki-ku, Yokohama, Kanagawa, Japan

